



Arms - Skin Cancer Excision



Forearm excision

Skin cancer on the forearm of an 85 year-old male following excision.



DermaClose application 2006 10

The DermaClose was applied in a radial configuration.



24 hours post application 2006

At 24 hours the wound had reduced in size approximately 90% and the tissue increased in laxity.



Linear Closure

The wound was sutured close 24 hours after application of DermaClose.

DERMAClose[®] RC
Continuous External Tissue Expander

WOUNDCARE[™]
TECHNOLOGIES
"Expanding Your Clinical Options"

Expanding your clinical options

www.dermaclose.com • (800)896-0436



Arms - Skin Cancer Excision cont.



At 8 weeks the wound has healed well.

 **DERMAClose[®] RC**
Continuous External Tissue Expander

WOUND CARE[™]
TECHNOLOGIES
"Expanding Your Clinical Options"

Expanding your clinical options

www.dermaclose.com • (800)896-0436